

SNELLEN TESTS

A few examples of Snellen tests available from Danker Lenses:




BC11931 Large Snellen "E" Test

20/200 to 20/10 size




BC 11942 *Tumbling E*

20 foot test card made of washable plastic.




BC1262 Illiterate E Chart



BCDI10 HOTV 10 Ft Test


This test for illiterate patients has HOTV

letters on one side and Allan figures on the other. Can be used with HOTV near card for matching with non-verbal patients.




BC814SCJ Snellen X-Cyl Jaeger Card

Near point test card printed on both sides. Front side tests for Jaeger 2 through 8. Back side tests Snellen and X Cylinder. Card is made of durable plastic and is washable.




BC1261 Snellen Test

With Red and Green lines to define 20/30 and 20/20 vision. 9" x 14"



BC11981 Near Point Test Card

Designed for use in performing the MEM (Monocular Estimate Method) test for determining lag of accommodation. The front of the card also had a cross-cylinder grid and a reduced Snellen chart. The back of the card provides additional reading material of various type sizes. The card is made of a plastic material which is easily cleaned.



BCNI HOTV Near Vision Acuity Card

This card is ideal for testing illiterate patients. It has HOTV letters on one side and Allen figures on the other