

DSA OPHTHALMIC STRIPS


Each strip is packed in paper pouch and sterilized by EO sterilization method.

DSA FLUORESCEIN:



- Standard (not high Molecular weight) for standard use and staining of the eye.
- Each strip contains Fluorescein Sodium U.S.P. 1 mg
- For Applanation Tonometry
- Can be used with RGP, Duette and Ultrahealth.
- NOT FOR NOVAKONE and SOFT LENSES
- Sterile individually wrapped Strips
- It is the most commonly used diagnostic strip used to visualize defects or abrasions in the corneal epithelium by staining the areas of cellular loss. Also useful for evaluating hard contact lens fitting & applanation tonometry.

Schirmer Tear Test Strips



This diagnostic strip made from the Whatman No 41 filter paper measures the amount of tear production. These strips are individually wrapped and are gauged so they can be easily read. The Schirmer test is based upon mechanical reflex stimulation of the tear production and is actually a test of reflex tear production when performed without anesthesia. Place a dry Schirmer strip over the middle to lateral 1/3 of the lower eyelid. The patient is then instructed to look forward and blink normally. The strip is removed after 5 minutes and amount of wetting is recorded in millimeters.

Product Details:

Model Name/Number	Schirmer Strips
Brand	Custom
Packaging Type	Box of 100 individual sterile strips
Type	Visual Acuity Examination Apparatus,Used in measuring amount of tear prodctcion
Number Of Strips	100 individual sterile strips
Material	Paper
Features	Used to Measure tear production

DSA Lissamine Green:

Lissamine Green Strips



These strips work by staining the cornea blue upon instillation. This causes "speckling" of the cornea, which gives the doctor immediate indication of any dry patches, as well as any muco-deficient or damaged corneal cells. Most importantly, it does not sting the eye like Rose Bengal, and a white light or blue light can be used on the slit lamp during detection. Available as sterile strips in individual pouches.

Box of 100 individually packed sterile strips.